


LimbLogic® Components

Vacuum Pump (LLV-21043*)



*Not sold individually. Available the following kits:

- LLV-2000-L LimbLogic Unilateral Kit, Lamination
- LLV-2002-L LimbLogic Bilateral Kit, Lamination
- LLV-2000-T LimbLogic Unilateral Kit, Thermoplastic
- LLV-2002-T LimbLogic Bilateral Kit, Thermoplastic

Charger (LLV-21011)




M6 X 55 mm Flat Head Capscrews, set of 4 (700-B092)



M6 X 50 mm Flat Head Capscrews, set of 4 (700-B093)

Exhaust Tube (LLV-21080)



Poron Filter (LLV-01071)



Proximal Sealing Gasket (LLV-01070)




Diagnostic Disc



Suction Pad (LLV-01100)



4-Hole Adapter (LLV-01041)



Vacuum Plate (LLV-01043)



Drop-In Adapter (LLV-02003)



Vacuum Pyramid (LLV-01044)

Optional:




4-Hole Adapter 1/8" Spacer Plate (700-300)

4-Hole Adapter 3/8" Spacer Plate (700-375)

Fob (LLV-21200)

Vacuum Pump Battery Indicator

Replacement Fob Battery:
commercially available
CR2032 coin cell battery



WillowWood®

The Ohio Willow Wood Company
15441 Scioto Darby Road
Mt. Sterling, OH 43143
phone 740.869.3377 / 800.848.4930
fax 740.869.4374 www.willowwoodco.com

Ohio Willow Wood Company B.V
Keizersgracht 62/64
1015 CS Amsterdam
The Netherlands



PN-2290-F 28 JAN 2014