

# ALPHA TRUSEAL LINER SIZE CHART

|                      | 25                               | 26.5                           | 28                              | 30                              | 32                              | 34                             | 37                               | 40                              | 43                             |
|----------------------|----------------------------------|--------------------------------|---------------------------------|---------------------------------|---------------------------------|--------------------------------|----------------------------------|---------------------------------|--------------------------------|
| Proximal Measurement | 37 - 49 cm<br>(14.5 - 19.25 in)  | 39 - 51 cm<br>(15.25 - 20 in)  | 40 - 54 cm<br>(15.5 - 21.25 in) | 42 - 58 cm<br>(16.5 - 23 in)    | 44 - 62 cm<br>(17.25 - 24.5 in) | 46 - 66 cm<br>(18 - 26 in)     | 49 - 72 cm<br>(19.25 - 28.25 in) | 52 - 78 cm<br>(20.5 - 30.75 in) | 55 - 85 cm<br>(21.5 - 33.5 in) |
| Distal Measurement   | 25 - 26.5 cm<br>(9.75 - 10.5 in) | 26.5 - 28 cm<br>(10.5 - 11 in) | 28 - 30 cm<br>(11 - 11.75 in)   | 30 - 32 cm<br>(11.75 - 12.5 in) | 32 - 34 cm<br>(12.5 - 13.5 in)  | 34 - 37 cm<br>(13.5 - 14.5 in) | 37 - 40 cm<br>(14.5 - 15.75 in)  | 40 - 43 cm<br>(15.75 - 17 in)   | 43 - 47 cm<br>(17 - 18.5 in)   |

  

The diagram shows nine Alpha TruSeal liner tubes of different sizes (25, 26.5, 28, 30, 32, 34, 37, 40, 43) arranged horizontally. Each tube is labeled with its size. Two horizontal dashed lines indicate measurement points: the proximal measurement is 20 cm (8 in) from the distal end, and the distal measurement is 4 cm (1.5 in) from the distal end. The tubes are shown in a perspective view, with the proximal measurement line at the top of the tube and the distal measurement line at the bottom of the tube.