

For order processing, ship order form and items to:  
 15441 Scioto Darby Road, Mt. Sterling, OH 43143 USA

### CLINIC INFORMATION

Clinic: \_\_\_\_\_

Prosthetist Name: \_\_\_\_\_

Email Address: \_\_\_\_\_ Date: \_\_\_\_\_

Bill to Account Number: \_\_\_\_\_ PO Number: \_\_\_\_\_

Shipping Address Line 1: \_\_\_\_\_

Shipping Address Line 2: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Preferred Contact Method: \_\_\_\_\_

Zip Code: \_\_\_\_\_ Phone: \_\_\_\_\_ Carrier:  UPS  FedEx

Shipping Method:  Ground  2-Day  Next Day

### PATIENT INFORMATION

Last, First Name: \_\_\_\_\_

Patient ID: \_\_\_\_\_ Amputation Level:  BK  Symes K-Level:  K1  K2  K3  K4

Left:  Right:  Bilat Order\*:

\*Bilat requires two order forms. Indicate side on each form

### DUPLICATE LINER INFORMATION

Duplicate Liner Serial Number \_\_\_\_\_ Qty \_\_\_\_\_

Cushion - ALC-DES-DP    Locking - ALC-DES-DP    DUO - ALC-DES-DP

### NEW LINER INFORMATION (Fill in for new liners or alterations to exiting liners)

MEASUREMENTS  centimeters  inches

**LINER TYPE**

Cushion  
ALC-DES-EO  
 Classic  
 Hybrid  
 Locking  
ALL-DES-EO  
 Classic  
 Hybrid  
 Duo  
ALD-DES-EO  
 Classic

**SHAPE CREATION**

Cast  
 STL / AOP  
 Omega / Tracer

CAD file name: \_\_\_\_\_

Cast Shipped Separately

**FABRIC TYPE**

Original  
 Buff  
 Gray/Green  
 Brown  
 Spirit  
 Buff  
 Gray/Green  
 Brown  
 Max - Buff  
 Select  
 Gray  
 Taupe

**LINER QUANTITY**

Qty \_\_\_\_\_

Additional Liner codes  
 Cushion    ALC-DES-DP  
 Locking    ALC-DES-DP  
 DUO        ALC-DES-DP

**LINER DESIGN**

Willow Wood to design liner based on shape and information provided (recommended)

Customer liner design based on notes and the thickness in each region indicated below

BASE LINER THICKNESS

mm    4mm to 12mm  
6mm Standard

DISTAL THICKNESS

mm    4mm to 25mm  
9mm Standard  
 4mm for Symes

PROXIMAL THICKNESS

mm    3mm to 9mm  
4mm Standard

POSTERIOR THICKNESS

mm    3mm to 9mm  
4mm Standard

FILL INVAGINATION

Yes (recommended)  
 No

STOVEPIPE LINER

Yes (bulbous distal)  
 No

Notes: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_