

Date MM / DD / YYYY

Please ship order form and items to:  
15441 Scioto Darby Road, Mt. Sterling, OH 43143 USA

**CUSTOMER INFORMATION**

Facility / Clinic: \_\_\_\_\_

Clinician Name: \_\_\_\_\_ Bill to Account #: \_\_\_\_\_ PO #: \_\_\_\_\_

Preferred Contact Method: \_\_\_\_\_ Shipping Address Line 1: \_\_\_\_\_

Email Address: \_\_\_\_\_ Shipping Address Line 2: \_\_\_\_\_

Phone: (\_\_\_\_) \_\_\_\_\_ - \_\_\_\_\_ City: \_\_\_\_\_ State: \_\_\_\_\_ Zip Code: \_\_\_\_\_

**SHIPPING INFORMATION**

Carrier:  UPS  FedEx      Method:  Ground  3-Day  2-Day  
 Next Day  Next Day Saver

**PATIENT INFORMATION**

Patient ID / Last Name: \_\_\_\_\_

Amputation Level:  TT  Symes      Amputation Side:  Left  Right  Bilateral

**\* Requires separate order forms. Indicate side on each form.**

**DUPLICATE LINER INFORMATION**

Duplicate Liner      Serial Number \_\_\_\_\_      Qty \_\_\_\_\_

Cushion - ALC-DES-DP      Locking - ALL-DES-DP      DUO - ADC-DES-DP

**NEW LINER INFORMATION** (Fill in for new liners or alterations to exiting liners)

MEASUREMENTS  cm  in

**LINER TYPE**

Cushion ALC-DES-EO  
 Classic  
 Hybrid  
 Locking ALL-DES-EO  
 Classic  
 Hybrid  
 Duo ADC-DES-EO  
 Classic

**SHAPE CREATION**

Cast  
 STL / AOP  
 Omega / Tracer

CAD file name: \_\_\_\_\_

Cast Shipped Separately

**FABRIC TYPE**

Original  
 Buff  
 Gray/Green  
 Brown  
 Spirit  
 Buff  
 Gray/Green  
 Brown  
 Max - Buff  
 Select  
 Gray  
 Taupe

**LINER QUANTITY**

Qty \_\_\_\_\_

Additional Liner codes  
Cushion ALC-DES-DP  
Locking ALL-DES-DP  
DUO ADC-DES-DP

**LINER DESIGN**

**WillowWood Design (Recommended)**  
Automatically based on cast or scan shape and order info. No contact required.

**Customer Design (Slower)**  
From notes & specific thicknesses indicated ↓

**BASE LINER THICKNESS**  
mm 4mm to 12mm  
6mm Standard

**DISTAL THICKNESS**  
mm 4mm to 25mm  
9mm Standard  
4mm for Symes

**PROXIMAL THICKNESS**  
mm 3mm to 9mm  
4mm Standard

**POSTERIOR THICKNESS**  
mm 3mm to 9mm  
4mm Standard

**FILL INVAGINATION**  
 Yes (recommended)  
 No

**STOVEPIPE LINER**  
 Yes (bulbous distal)  
 No

Notes: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_